

FLUORIDE PROBES & SUPPLIES

THERMO SCIENTIFIC ORION FLUORIDE IONPLUS COMBINATION PROBE

For rapid determinations of Fluoride in water. **Polyphosphates do not interfere.** No separate reference probe required. Includes probe and one 60 ml bottle of filling solution. BNC connector.

O-9609BN 765.00/ea

FILLING SOLUTION FOR O-9609BN ONLY

O-900061 21.40/ea 99.00/pk of 5 x 60 ml



THERMO SCIENTIFIC ORION FLUORIDE PROBE. BNC connector. Requires O-900100 Reference probe.

Description	NCL #	Price/ea
Fluoride probe	O-9409BN	657.90
Single junction reference probe	O-900100	217.55

Filling Solution for O-900100 when used with O-9409BN Fluoride Probe **ONLY**
O-900001 19.75/ea 90.00/pk of 5 x 60 ml

ISA Buffers for use with Fluoride Probes

Please see the chemical section in this catalog for Fluoride buffers in other sizes.

Concentration	Manufactured by	NCL #	Price
TISAB III Concentrate	Thermo Scientific Orion	O-940911	145.35/475 ml HM
TISAB III Concentrate	NCL	F-29	37.35/500 ml HM
TISAB II	Thermo Scientific Orion	O-940909	136.80/3.8 liter
TISAB II	NCL	F-29A	75.60/3.8 liter

Note: TISAB II is mixed 1:1 with samples and standards. TISAB III is mixed 1:10 with samples and standards. Therefore, TISAB III would last you 10 times as long as TISAB II.

Fluoride Standards

Please see the chemical section in this catalog for Fluoride standards in other sizes.

Concentration	Manufactured by	NCL #	Price
100 ppm as F	Thermo Scientific Orion	O-940907	71.60/475 ml
100 ppm as F	NCL	F-39	20.60/500 ml
10 ppm as F with TISAB II	Thermo Scientific Orion	O-040908	84.00/475 ml
10 ppm as F with TISAB II	NCL	F-39Q	27.45/500 ml
2 ppm as F with TISAB II	Thermo Scientific Orion	O-040907	84.00/475 ml
2 ppm as F with TISAB II	NCL	F-39P	27.45/500 ml
1 ppm as F with TISAB II	Thermo Scientific Orion	O-040906	84.00/475 ml
1 ppm as F with TISAB II	NCL	F-39M	27.45/500 ml